

Checklist

Preventive Maintenance for A-dec 41L/42L Simulators

A-dec recommends that you perform preventive maintenance annually or after 1500 hours of use.

General

Notes

- Check the operation of the torso movement including up, down, and rotation.
- Use the touchpad (if equipped) to check for proper function and advise on recommended repairs or replacements.

Delivery System

Notes

- Use the touchpad (if equipped) and the foot control to check the operation of handpiece identification and activation, if applicable.
- Use the syringe nut torque adaptor (p/n 23.1354.00) to torque the syringe nut.
- Check for damaged, brittle, or stiff handpiece tubing and recommend replacement as necessary.
- Check the condition of the handpiece coupler O-rings and advise on repair as necessary.
- Check the control block valve condition. Replace the valve stem O-rings every five years or more frequently, if needed.
- Check the handpiece pressures and adjust to manufacturer's specification as necessary.
- Check the fiber optic voltages (if equipped) and adjust to handpiece manufacturer's specification as necessary.
- Replace the syringe O-rings and button-valve assemblies, lubricate with A-dec silicone lubricant, and verify the operation.
- Check for oil or water accumulation in the oil collector, replace gauze, and clean as necessary.
- Replace the self-contained water bottle sealing washer or O-rings (Self-Contained Water Option Only).
- Test the function of all touchpad buttons.
- Check the operation of the handpiece holder valves and repair/replace as required.

Foot Control

Notes

- Check the condition of foot control tubing.
- Inspect and replace O-rings and the trim ring as necessary. Lubricate O-rings with A-dec silicone lubricant, and check for correct operation.
- Replace diaphragm every five years or more frequently, if needed.

Utilities

Notes

- Replace air filter (p/n 24.0234.00) and water filter (p/n 41.0004.00) screens.
- Inspect and replace air and water regulator diaphragms and O-rings as necessary.
- Check incoming air and water pressures and adjust as necessary.
- Check for air and water leaks.
- Check for damaged, brittle, or stiff tubing.
- Check incoming air to ensure clean and dry air is being delivered to the delivery system.
- Check moisture separators at the compressor, and replace moisture separator filters as necessary.

Other

Notes

- Check the monitor mount tilt tension and adjust as necessary.
- Check all arm pivot points for correct tension and adjust as necessary.
- Inspect and adjust fasteners as needed.
- Inspect covers to ensure proper fit and adjust as needed.
- Check caster rotation to ensure debris is not inhibiting movement.
- Check stool cylinder for leakage and ensure rod is free of visible dirt/grease.
- Check all stool seat, back/torso, and height adjustments to ensure proper operation.
- Check upholstery condition and recommend replacement as appropriate.
- Touch up “through to metal” scratches with paint as necessary.

Review

Notes

- Ensure staff and students are aware of the location of all Instructions for Use documents for equipment.
- Report on any necessary repairs and provide quotes as required.
- If required, instruct staff on correct maintenance and operation.
- Review A-dec’s Smiluator *Equipment Maintenance Checklist* (p/n 86.0703.00) with staff to ensure that routine maintenance is occurring as required. Point out checklist items that you have just completed as part of this preventive maintenance.

School / Clinic: _____
Address: _____

Station #: _____

Service Technician: _____
Service Date: _____

Reviewed with staff / customer:

Signature and Date

Staff / customer unavailable to review checklist

Additional Comments / Recommendations: _____



A-dec Customer Headquarters

Tel: 1.800.547.1883 within USA/CAN
Tel: +1.503.538.7478 outside USA/CAN

A-dec China

Tel: +1.503.538.7478

A-dec Australia

Tel: 1.800.225.010 within AUS
Tel: +61.(0).2.8332.4000 outside AUS

A-dec United Kingdom

Tel: 0800.ADEC.UK (2332.85) within UK
Tel: +44.(0).24.7635.0901 outside UK



86.0704.00 Rev A
Copyright 2017 A-dec Inc.
All rights reserved.