

Cavitron.

For dental professionals searching for the streamlined look of integration, **Cavitron® Integrated Systems** gives you the freedom of more open space.

Cavitron® Integrated Systems

- Fewer cords and less visible equipment
- Open office environment preserves counter space
- Compatible with Steri-Mate® 360° Rotating Handpiece



There's only one Cavitron®.

Cavitron.
integrated
systems

Cavitron® Integrated Systems offer an open environment with the same benefits of traditional Cavitron® units.

- Compatible with **Steri-Mate® 360°**, a rotating handpiece that allows the clinician to fully rotate an ultrasonic insert without stopping for adjustments
- Detachable and autoclavable handpiece reduces cross contamination for enhanced infection control
- **Patented Sustained Performance System, SPS™, Technology** automatically adjusts power to maintain scaling efficiency, even at low power
- **BlueZone™** extended low power range designed for improved patient comfort



Cavitron Integrated Systems integrate easily with A-dec® 500 and A-dec® 300 delivery systems.

For dental professionals searching for the streamlined look of integration, you are now able to combine Cavitron Systems with the convenience of integration into the A-dec 500 or A-dec 300 dental units. You can also put more technology at your fingertips with the A-dec deluxe touch pad, which gives you an intuitive display, easy-touch control, and programmable memory for preferred ultrasonic power settings.

Item	A-dec Integrated Part #	Description	Compatible With:
Cavitron A-dec Built-In Ultrasonic Scaler Kit	54.0599.00	Cavitron Steri-Mate 360° Kit, 61"	A-dec 532, 533, 542, 332, 333, 334, 335, 336
	54.0597.00	Cavitron Steri-Mate Kit, 61"	
	54.0600.00	Cavitron Steri-Mate 360° Kit, 82"	A-dec 541
	54.0598.00	Cavitron Steri-Mate Kit, 82"	

For an integrated Cavitron kit contact:

1.800.547.1883 Within USA
1.503.538.7478 Outside USA/Canada
www.a-dec.com



For a built-in Cavitron kit for field installation contact:

DENTSPLY Professional
1.800.989.8826
www.dentsply.com

