

A-dec® 500 Fleet Packages

For Dental Service Organizations

- Chair
- Delivery
- Stools
- Assistant's Package
- LED Light



Assistant's package not shown

Contact the A-dec DSO team today. Call 1.800.547.1883.

Make a smart investment in your future with A-dec solutions that help you work more productively, comfortably, and intuitively. Configurable for your needs, A-dec 500 Fleet packages start with our most comfortable chair, and offer a range of customizable options—together with legendary A-dec quality and reliability.

Every A-dec 500 Fleet package includes:

- Dental chair (511ft)
- Traditional delivery system (332ft)
- Doctor's and assistant's stools (521 and 522, or 421 and 422)
- Radius® LED light (372L or 572L)
- Assistant's package (351)

511ft chair features:

- Pressure-mapped surfaces and 9" toeboard tilt—for best-in-class patient comfort
- Virtual pivot—synchronizes with natural movement of the patient
- Low base-down point and ultra-thin, flexible backrest—for optimal access
- 60° chair swivel with brake—for easy rotation and patient positioning
- Dual-articulating, gliding headrest—comfortably supports patient's head and neck
- Precision hydraulic lift and tilt system—for smooth, quiet, quick positioning
- Contoured floor box—discreetly houses air and water controls
- 300-watt power supply

332ft delivery system features:

- Balanced flex arm with air brake—for easy positioning and stability
- Four handpiece holder positions—customizable for multiple tools
- Radius-mounted—rotates around the dental chair for right- and left-hand compatibility
- 2-liter, quick disconnect, self-contained water bottle—to easily control treatment water quality
- Adjustable air coolant—prevents handpiece damage from frictional heat
- Individually adjustable water coolant—accommodates coolant needs for each handpiece

351 Assistant's package features:

- Autoclavable HVE
- Saliva ejector
- Syringe


Customize your 500 Fleet package



Sewn or formed upholstery



400 or 500 level stools



3 or 4 position control block



372L or 572L dental light



Factory-integrated clinical devices, including electric motors and scalers



Standard or Deluxe touchpad



a-dec.com/DSO-Solutions
DSOcustomerservice@a-dec.com

The colors and products shown may vary from the actual materials and are subject to change without notice. Please see your authorized A-dec dealer for accurate color samples and the most current product information.