

# Intego

A-dec Signature  
Handpieces and Motors



# Reliable performance, impressive value.

You know exactly what you need for your practice. That's the inspiration behind our signature line of A-dec handpieces and motors.

Enjoy the thoughtful design and quality craftsmanship you expect, backed by legendary customer service and support. This singular line of motors and attachments is an ideal solution for your essential, day-to-day procedures. And an easy choice to make when you're selecting your A-dec chair, lights, and delivery.

You earn your patients' trust with every treatment. And we'll earn yours every time you reach for an A-dec handpiece.



Not available in all regions. Check with your authorized A-dec dealer for local availability.

# CLASS A A-DEC LOW-SPEED AIR COLLECTION

A-dec low-speed air handpieces offer durability and reliable performance. Their sleek, lightweight design ensures an entirely comfortable balance in your hand. High-precision gears provide quiet operation and reduce uncomfortable vibration. The extra-slim head and body improve visibility as you work, even in difficult-to-access areas of the oral cavity.

## Designed for long life.

The debris defense system on A-dec low-speed air handpieces automatically reduces the entry of oral fluids and contaminants into the handpiece head to prolong the life of the bearings. The prophy handpiece features a paste seal system to reduce dust and polishing paste debris from entering the handpiece; and a special rotor coating prevents potential breakage caused by surface abrasions from the polishing paste.



### A11

<b>Application</b>	General application handpiece
<b>Gear Ratio</b>	1:1 direct drive
<b>Max Speed</b>	40,000 rpm
<b>For</b>	Bur for Contra Angle (CA) 2.35mm diameter
<b>Head Size</b>	Φ9.6mm×13.5mm
<b>Warranty</b>	1 year

### Features

- Non-optic
- Debris defense system
- Push button chuck
- Ball bearing
- Aluminum body
- Autoclavable



### A41

<b>Application</b>	General application handpiece
<b>Gear Ratio</b>	4:1 reduction
<b>Max Speed</b>	10,000 rpm
<b>For</b>	Bur for Contra Angle (CA) 2.35mm diameter
<b>Head Size</b>	Φ9.6mm×13.5mm
<b>Warranty</b>	1 year

### Features

- Non-optic
- Debris defense system
- Push button chuck
- Ball bearing
- Aluminum body
- Autoclavable



### A11S

<b>Application</b>	General application handpiece
<b>Gear Ratio</b>	1:1 direct drive
<b>Max Speed</b>	40,000 rpm
<b>For</b>	Bur for Handpiece (HP) 2.35mm diameter
<b>Head Size</b>	N/A
<b>Warranty</b>	1 year

### Features

- Non-optic
- Debris defense system
- Aluminum body
- Autoclavable



### A41P

<b>Application</b>	Prophylaxis handpiece
<b>Gear Ratio</b>	4:1 reduction
<b>Max Speed</b>	5,000 rpm
<b>For</b>	Screw in cups and brushes
<b>Head Size</b>	Φ7.4mm×12mm
<b>Warranty</b>	1 year

### Features

- Non-optic
- Paste seal system
- Aluminum Body
- Autoclavable

## A-DEC 301 AIR MOTOR

With its clean aluminum design, the A-dec 301 motor offers a balanced combination of light weight and durability. Enjoy smooth forward and reverse drive functionality, with an adjustable speed range up to 22,000 rpm. The A-dec 301 motor accommodates non-optic handpieces and contra-angles with E-type compatibility.

<b>Optic</b>	Non-optic
<b>Speed</b>	22,000 rpm +/- 10%
<b>Weight</b>	62 grams
<b>Material</b>	Aluminum
<b>Tubing</b>	Midwest 4-hole
<b>Warranty</b>	1 year
<b>Autoclavable</b>	Yes
<b>Water Spray</b>	No



Not available in all regions. Check with your authorized A-dec dealer for local availability.

## A-DEC 501 ELECTRIC MOTOR

Responsive and lightweight, the A-dec 501 electric motor delivers all the power you need, with a full range of speed. Work with more precision and consistency—and less fatigue. Program the settings to match your preferences for activation, speed, torque and direction.

<b>Applications</b>	General restorative, prophylaxis, endotics	<b>Weight</b>	57 grams (micromotor only)
<b>Contra-angle compatibility</b>	E-type	<b>Optic</b>	Daylight quality LED
<b>Motor Type</b>	Brushless micromotor	<b>Warranty</b>	1 year
<b>Speed (Standard Mode)</b>	1,000 - 40,000 rpm	<b>Autoclavable</b>	Yes
<b>Speed (Endo Mode)</b>	100 - 5,000 rpm	Couple the A-dec 501 electric motor with the A-dec Deluxe touchpad to fully integrate your delivery system—and add the capability to perform endodontic procedures and functions such as torque control and automatic reverse.	
<b>Torque (Standard Mode)</b>	Max. 4.0 Ncm	An electric motor controller is a required component of the A-dec delivery system, which is configured at the time of your A-dec purchase.	
<b>Torque (Endo Mode)</b>	Max. 3.0 Ncm		



Not available in all regions. Check with your authorized A-dec dealer for local availability.

**A-dec Headquarters**

2601 Crestview Drive

Newberg, Oregon 97132 USA

T: 1.800.547.1883 within USA/CAN

T: +1.503.538.7478 outside USA/CAN

F: +1.503.538.0276

**a-dec.com**

**A-dec Australia**

T: 1.800.225.010 within AUS

T: +61.(0)2.8332.4000 outside AUS

**A-dec China**

T: 400.600.5434 within CN

T: +86.571.89026088 outside CN

**A-dec United Kingdom**

T: 0800.233285 within UK

T: +44.(0).24.7635.0901 outside UK



The products shown may vary and are subject to change without notice. Contact your authorized A-dec dealer for the most current product information. All trademarks in this document are subject to the trademark rights of A-dec, its affiliates, and their collective suppliers in the United States and foreign countries.

©2021 A-dec Inc. All rights reserved.

85.1021.00/BB/3M/7-21/(9-22)/Rev A2