

Convenience and Flexibility

Cascade Carts



 **a dec**[®]

Flexible, stable, and reliable delivery.



Cascade duo cart



Cascade assistant's cart



Cascade cart

Flexible and convenient choices

A-dec's Cascade® carts offer you an excellent combination of choices, flexibility, and convenience to streamline your operatory and increase efficiency. You can choose three different styles: the Cascade Duo Cart, which features doctor's and assistant's instrumentation; the Cascade Assistant's Cart with integrated assistant's instrumentation; and the Cascade Cart. All three carts feature a sturdy, height-adjustable worksurface that (on duo and assistant's carts) offers you a choice from a full range of laminate and solid surface colors.

For flexible and convenient instrumentation delivery, choose the duo or assistant's carts. In addition to the worksurface, both carts include a rotating arm that delivers assistant's instrumentation precisely where you need it. On the duo cart, the control unit rotates as well, so both carts can accommodate either right- or left-handed operators. The utilities of these carts can be mounted to the wall or in your cabinetry, allowing you more room to move around. And with four durable, smooth-rolling casters, you can easily position your cart during procedures and quickly move it out of the way when it's not in use.

Reliability you can count on

Cascade carts are built to last. Duo carts, for instance, feature A-dec's proven, reliable Century Plus controls with doctor-replaceable cartridges, so you can easily keep your system operating like new. And since all A-dec carts feature durable compounds, stain-resistant materials, and rugged components such as a solid steel base to maximize stability, you can be sure your Cascade cart will easily withstand the daily rigors of your practice for years to come.

Standard Product Features	Cascade Duo	Cascade Asst's
Pivoting Cascade delivery system with Century Plus control block	●	
Handpiece oil collection system	●	
Wet/dry foot control and carrier	●	
Cascade solids collector	●	●
Pivoting instrumentation arm with autoclavable HVE and saliva ejector	●	●
Self-contained water system	●	●
Water QD outlet and flow control	●	●
Cable prewiring	●	●
Duplex electrical outlet	●	●
Height-adjustable worksurface	●	●
Laminate worksurface	●	●

Optional Product Features	Cascade Duo	Cascade Asst's
Autoclavable syringe (<i>doctor-side</i>)	○	
Autoclavable syringe (<i>assistant's-side</i>)	○	○
QDs for air or self-contained water	○	○
Dual HVE	○	○
Single or dual voltage intraoral light source	○	
Programmable chair touchpad	○	○
Arm-mounted tray holder	○	○
Solid surface or chemical resistant laminate worksurface	○	○

A-dec Inc.
2601 Crestview Drive
Newberg, OR 97132 USA
Phone: 1.800.547.1883, 1.503.538.7478
Fax: 1.503.538.0276
www.a-dec.com

A-dec, Cascade, Decade, Performer, and the A-dec logo are registered trademarks in the U.S. Patent and Trademark office. All other products or services mentioned in this document are covered by the trademarks, service marks, or product names as designated by the companies that market those products. Touchpad symbols are proprietary to A-dec Inc.

The colors and products shown may vary from the actual materials and are subject to change without notice. Please see your authorized A-dec dealer for accurate color samples and the most current product information.