

Setting the Standard

Cascade Chair



 **a dec**[®]

Quality, style and reliability



Thin-profile backrest casting



Automatic armrests



Durable iron and aluminum castings

Comfort for you and the patient

The Cascade[®] 1040 chair is designed to offer a balance of both comfort and access to the patient. Cascade's thin-profile backrest casting, for example, allows you to get your legs comfortably under the chair. To keep thickness to a minimum while offering comfort, the chair upholstery features layers of different foam densities—one for lumbar support, one for overall support, and one for contoured comfort. And to provide even more cushioning for your patients, you can choose to have your Cascade chair fitted with optional sewn upholstery.

Rugged and Reliable

Look beneath the surface of the chair, and you'll begin to discover even more evidence of Cascade's tremendous value. Like ten rugged castings that ensure a stable structure, a simple microprocessor that allows you to program the chair without complication, and an ultra-quiet hydraulic system with proven reliability. So whether you're looking at Cascade's streamlined exterior, or examining Cascade's durable construction, you can be sure every detail will maximize your investment and help you provide your patients with the finest in dental health care.

For a closer look at the Cascade chair, contact your local authorized A-dec dealer today, or call us at 1.800.547.1883.

Product Features	Standard	Optional
Sewn Upholstery		<input type="radio"/>
Programmable button footswitch	<input checked="" type="radio"/>	
Programmable chair touchpad (<i>on delivery system</i>)		<input type="radio"/>
Thin-profile backrest casting	<input checked="" type="radio"/>	
Lumbar support	<input checked="" type="radio"/>	
Automatic armrests	<input checked="" type="radio"/>	
Cast-aluminum baseplate	<input checked="" type="radio"/>	
Double-articulating headrest	<input checked="" type="radio"/>	
Hydraulic-lift cantilever system	<input checked="" type="radio"/>	
Chair swivel with mechanical brake	<input checked="" type="radio"/>	
Duplex electrical outlet	<input checked="" type="radio"/>	
Internal umbilical for integrated Cascade or Radius [®] delivery system		<input type="radio"/>
Chair mounts for Radius delivery system, flat-panel monitor mount, and dental light		<input type="radio"/>