

Radius[®] on A-dec 511



Combine the ever-popular Radius delivery system with the award-winning A-dec 511 chair, and you have an exceptional, "entry-level" option for integrating technology into your treatment rooms.


Another Quality Choice From A-dec



Features and Options for Radius on A-dec 511



Radius Traditional Delivery



Radius Continental® Delivery



A-dec 561 Cuspidor

When combining the award-winning A-dec 511 chair with Radius delivery, you get an incredibly comfortable chair surrounded by a full range of quality of choices, including two styles of delivery, cuspidor, multiple configurations of assistant's instrumentation, dental light, monitor mounting system, and many more options.



A-dec 551 Long Assistant's Arm



A-dec 551 Short Assistant's Arm



A-dec 571 Dental Light



A-dec 561 Side-Support Monitor Mount

- **Seamless or Sewn Upholstery** on the A-dec 511 chair with built-in lumbar support and memory foam to enhance patient comfort
- **Rotating Delivery Arm** offers excellent positioning and full left/right compatibility
- **Self-Contained Water System**, combined with **A-dec ICX® Waterline Tablets**, for greatly simplified waterline maintenance
- **AlphaSan® Antimicrobial Tubing** throughout the system and handpiece tubings to further protect against waterline contamination
- **Rotating Assistant's Arm or Cuspidor** for left/right compatibility, exceptional positioning, and improved ergonomics
- **Dental Light or Monitor Mount Arms** pivot around the back of the chair to accommodate left- and right-handed operators
- **Assistant's Instrumentation** includes fully autoclavable syringe, HVE, and saliva ejector to help simplify asepsis
- **Integrated Assistant-Side Touchpad** serves as a convenient secondary control for basic chair, light, and cuspidor functions
- **Optional A-dec Intraoral Camera** with live, freeze, and quad split-screen views activated through Radius foot control
- **Optional A-dec Electric Micromotor** offers brushless technology for less maintenance and a much longer life cycle than previous micromotor technology



A-dec Inc.
2601 Crestview Drive
Newberg, Oregon 97132 USA
Phone: 1.800.547.1883, Fax: 1.503.538.0276
Outside USA/Canada: 1.503.538.7478

For more information about Radius on A-dec 511, contact your local authorized A-dec dealer, visit www.a-dec.com, or call us at 1.800.547.1883 (1.503.538.7478 worldwide).

The colors and products in this brochure may vary from the actual materials and are subject to change without notice. Please see your authorized A-dec dealer for accurate color samples and the most current product information.

© 2006 A-dec® Inc. All rights reserved.
85.6113.00 / DYN / 10M / 02-06

www.a-dec.com