



The A-dec 500® deluxe touchpad offers manual and automatic settings, including up to 16 programmable memory presets. (8 in standard mode and 8 in endo mode)

the strong, versatile type

It's one of the smallest, lightest electric motors ever made. The ultra-precise EA-51LT, brought to you in partnership by A-dec and W&H, delivers unprecedented power, with unparalleled precision and control. Fully integrated into the A-dec 500, the EA-51LT provides over 60 watts of high-torque power to your chosen attachment in two operative modes. Whether for standard restorative to sophisticated endodontic work, now you can take it all on with total confidence.

Smarter for your practice, safer for your patients

The EA-51LT offers two operating modes and a long list of features:

- Standard and endodontic alternated by the push of a button
- Over 60 watts of high-torque power
- Ergonomic design
- Touchpad and/or foot control for uninterrupted focus
- Ultra quiet operation
- Seamless integration with the A-dec 500 eliminates downtime
- 3 to 200,000 rpm with A-dec/W&H attachments
- Torque control with auto-reverse



Maintains high-torque power at speeds between 3 and 200,000 rpm for greater control and accuracy in all applications

the motor behind the handpiece

With A-dec/W&H attachments, the EA-51LT can maintain high-torque power at speeds as low as 3 rpm and as high as 200,000 rpm. That means greater control and accuracy in your hands – from high-speed cutting to refining and endodontics. So now you'll have the control and the power you need for perfect procedure after perfect procedure.

Legendary reliability

The EA-51LT was built for A-dec by W&H, a company with a reputation for creating the finest, most trusted dental products on the planet. W&H is one of the world's largest producers of rotary instruments and has been a family-owned business for more than 100 years. Together, A-dec and W&H offer the newest evolution of brushless technology. When you choose the EA-51LT, you can expect incomparable reliability and reduced maintenance concerns to keep your practice productive and profitable.



A-dec Inc.
2601 Crestview Drive
Newberg, Oregon 97132 USA
Phone: 1.800.547.1883, Fax: 1.503.538.0276
Outside USA/Canada: 1.503.538.7478

© 2006 A-dec® Inc. All rights reserved.
85.6115.00 / ??? / 40M / 04-06

To find the full list of benefits of the EA-51LT and A-dec's versatile endo package, contact your local authorized dealer, call 1.800.547.1883 or visit www.a-dec.com.

www.a-dec.com