

Get to know our family

A-dec's complete line of healthy dental equipment solutions creates a collection of possibilities for a treatment room that works exactly how you want it to. Chair, delivery and support systems combine to give you unparalleled range of motion, industry-leading technology accommodation, plus greatly simplified maintenance and asepsis innovations. These are just a few examples of how A-dec has redefined quality dental equipment — and provides the best return on your investment.



World Headquarters	International Distribution Centers	
2601 Crestview Drive	A-dec Australia	A-dec United Kingdom
Newberg, Oregon 97132 USA	Tel: 1.800.225.010 Within Australia	Tel: 0800 ADECUK (233285) Within UK
Tel: 1.800.547.1883 Within USA/Canada	Tel: + 61 (0)2 8332 4000 Outside Australia	Tel: +44 (0) 24 7635 0901 Outside UK
1.503.538.7478 Outside USA/Canada	a-dec.com.au	a-dec.co.uk
Fax: 1.503.538.0276		
a-dec.com		



A-dec, A-dec logo, A-dec 300, A-dec 500, Continental, Decade, Radius, ICV, Preference, Preference Collection, and Preference ICC are registered trademarks in the U.S. Patent and Trademark Office. All other products or services mentioned in this document are covered by the trademarks, service marks, or product names as designated by the companies that market those products. Touchpad symbols are proprietary to A-dec Inc. The colors and products shown may vary from the actual materials and are subject to change without notice. Please see your authorized A-dec dealer for accurate color samples and the most current product information.
©2010 A-dec Inc. All rights reserved. 85.6157.00/GAC/19M/05-10/Rev B



premium

A-dec 500[®]




classic

Cascade[®] 1040



Product Solutions



simple

Decade[®] Plus



compact

A-dec 300[®]



A-dec 500

Advanced innovation. A-dec 500 delivers optimal performance and patient comfort for the doctor who understands the value of progress. With outstanding access to the patient, seamless integration of accessories, and anatomical contours to ensure patient comfort, this is the hallmark of treatment room excellence. With its streamlined profile and complete range of features, A-dec 500 is unmatched in its elegant display of innovation.

A-dec 500 – Great ergonomic access for doctor, comfort for patient and full internal integration of all ancillaries.

FEATURE	BENEFIT
Ultra-thin 1" (25 mm) backrest	Optimal oral cavity access Increased leg room under the chair
Seamless ancillary integration	Allows equipment to meet the changing needs of dentistry
Anatomically shaped contours	Reduces pressure points Comfortable and supportive
Unique self-contained water bottle	Reduces risk of cross-contamination
Exceptional range of motion plus low base-down position (13.5"/343 mm)	Improves access, visibility and posture
Integrated touchpad	Easy control of chair, cuspidor, dental light
Rotating arm system and pivoting, centrally mounted control head	Advanced instrument positioning Extended tubing reach

Delivery & Options

- A-dec 532 Traditional
- A-dec 533 Continental
- A-dec 551 assistant's instrumentation
- A-dec 561 support center, cuspidor & monitor mount
- A-dec 571 dental light
- Upholstery: seamless or sewn



Cascade 1040

Cascade is the solution that has a way of putting the patient first without compromising the ease of the dentist. Layers of various foam densities – one for lumbar support, one for overall support, and one for contoured luxury – ensure long-lasting patient comfort, while easy-to-program presets and reliable ultra-quiet hydraulics make Cascade 1040 a satisfying plus for any treatment room. (Note: Intraoral cameras, electric motors, curing lights, and ultrasonic handpieces require an external mount.)

Cascade 1040 – Patient comfort and positioning for long procedures and delivery system with external integration of ancillaries.

FEATURE	BENEFIT
Sculpted backrest	Patient comfort
Armrests tuck into chair	Easy patient entry/exit
Smooth rounded surfaces	Quick and simple asepsis
Double-articulating headrest	Adjusts patient for better access and viewing of oral cavity
Programmable button footswitch or optional touchpad	Allows manual or automatic chair positioning
Optional internal umbilical and integrated Radius delivery system	Frees up treatment room space

Delivery & Options

- Radius 2122 Traditional
- Radius 2132 Continental
- Radius 7115 assistant's instrumentation with integrated touchpad
- Radius 7285 cuspidor
- Radius dental light
- Upholstery: seamless or sewn



Decade 1221

Uncomplicated ease. With such a high level of proven reliability and range of options you'd expect on higher-priced equipment, Decade Plus 1221 delivers simple-to-operate versatility. And with its double-articulating headrest, quick-release armrests and 60° swivel, there's no sacrificing easy patient accessibility or comfort. Designed to maximize productivity, this is the solution to work with your existing delivery system or integrate with the left/right compatibility of Radius® for one of the best chair values available.

Decade Plus 1221 – fully positionable chair and delivery system with external integration of ancillaries.

FEATURE	BENEFIT
Thin 1-1/2" (38 mm) backrest	Easy patient access Comfortable working position
Soft cushioned seat, smooth motion	Enhanced patient comfort
Supportive quick-release armrests	Easy entry/exit and stand-up dentistry
Double-articulating headrest	Adjusts patient for better access and viewing of oral cavity
Compatible with chair adapters	Accommodates existing delivery system
60° swivel (30° each side)	Works behind the patient even in small treatment rooms
Optional internal umbilical and integrated Radius delivery system	Frees up treatment room space

Delivery & Options

- Radius 2122 Traditional
- Radius 2132 Continental
- Radius 7115 assistant's instrumentation with integrated touchpad
- Radius 7285 cuspidor
- Radius dental light
- Upholstery: seamless or sewn



A-dec 300

Simple. Smart. Streamlined. A-dec 300 takes innovative design to a whole new level. Stylish yet simple. Although A-dec 300 may be compact in size, it's big in versatility and integration. By combining the Radius®-style delivery options and left/right assistant's arm – including a wide range of ancillary possibilities – this is the smart, affordable system that adapts to whatever your future brings. As you grow, so can your dental equipment.

A-dec 300 – Sleek, stylish design available in seamless or sewn upholstery and with a complete range of integrated options.

FEATURE	BENEFIT
Ultra-thin 1-1/2" (38 mm) backrest	Good oral cavity access, improved posture
Radius-style delivery option	Versatile left/right compatibility
Chair-base mount support center	Small, sleek footprint
13.75"-29.5" (349 - 749 mm) vertical moving range	Accommodates stand-up dentistry
Simplified design with fewer moving parts	Lighter weight Simplifies maintenance and cleaning Increases durability
Open toeboard	Cleanable
Patient-adjustable neck support option	Lets the patient self-adjust the neck-support
Double-articulating headrest option	Adjusts the patient for better access and viewing of the oral cavity
Seamless ancillary integration	Allows equipment to meet the changing needs of dentistry

Delivery & Options

- A-dec 334 base or post-mount Traditional
- A-dec 332 Radius-style Traditional
- A-dec 335 base or post-mount Continental
- A-dec 333 Radius-style Continental
- A-dec 351, 352 or 353 assistant's instrumentation
- A-dec 361 cuspidor and support center
- A-dec 371 dental light
- Upholstery: seamless or sewn