

EA-41

Electric Motor



COMPATIBLE. Versatile tabletop configuration works with your existing delivery system

COMPACT. Smallest, lightest electric motor available

INNOVATIVE. Superior LED lighting available with A-dec |W&H Alegra attachments

AFFORDABLE. All at a price that's in your budget



The optimal choice for restorative dentistry

Get the size and shape of an air-driven highspeed *and* the power and precision cut of an electric highspeed or lowspeed

EA-41 Specifications:

	EA-41
Type of motor	Electric
Weight	2.25 oz (64 g)
Length	38.55 mm
Attachment coupling	E-style
Light	No (connects to A-dec W&H Alegra attachments that have their own self-generated LED light)
Speed range (RPM)	2,000-40,000
Maximum power	55 W
Torque	3 Ncm
Delivery system	A-dec Radius or tabletop
Rotary direction	Forward and reverse



Features:

- 20% shorter than any other motor (38.55 mm)
- 16% shorter length than typical electrics
- Lightest electric motor available (64 g)
- LED lighting from Alegra attachments
- Highspeed and lowspeed applications
- 2,000-40,000 RPM in forward and reverse directions for general restorative functions
- E-coupling (use with A-dec|W&H Alegra attachments with LED light)
- Connects to any other E-style attachment (will not provide light)
- Sterilizable
- Third-generation proven brushless motor technology

Benefits:

- Versatile—designed to work with your existing delivery system or A-dec Radius® delivery systems.
- Compact and uniquely balanced to reduce hand strain and increase maneuverability.
- The EA-41 motor is available with eight attachment options (highspeed and lowspeed contra-angles, straight).



Get LED lighting when you use Alegra attachments

The W&H Alegra Handpieces shown in this brochure are available through A-dec only in the USA, Canada, and Australia. For all other countries, please visit www.wh.com to find your local dealer.

A-dec Headquarters
2601 Crestview Drive
Newberg, Oregon 97132 USA
Tel: 1.800.547.1883 Within USA/Canada
1.503.538.7478 Outside USA/Canada
Fax: 1.503.538.0276
www.a-dec.com

International Distribution Centers

A-dec Australia
Tel: 1.800.225.010 Within Australia
Tel: + 61 (0)2 8332 4000 Outside Australia
www.a-dec.com.au

A-dec United Kingdom
Tel: 0800 ADECUK (233285) Within UK
Tel: +44 (0) 24 7635 0901 Outside UK
www.a-dec.co.uk

