

# A-dec 300 Package

For Dental Service Organizations



 a dec<sup>®</sup>

We've distilled the most essential features into one smart solution. Compact. Sensible. Flexible. A-dec 300 brings you exactly what you need with an added bonus—it's from A-dec, so you know it's built to last.

## Features and Benefits

- 60° chair swivel with brake—allows easy rotation and patient positioning.
- Dual-articulating, gliding headrest—supports head and neck, moves with the patient, and adjusts with one hand.
- Fully ambidextrous delivery system—change from right-hand to left-hand configuration.
- Hydraulic lift system—A-dec built for robustness and gentle stops and starts during patient positioning.
- Large 2-liter water bottle—requires fewer refills.
- Safety lift arm stop plate—stops immediately on contact for operator safety.
- Contoured floor box—sleek, unobtrusive floor box discreetly houses air and water controls.

## A-dec 300 DSO Package Includes

- Chair with Contoured Floor Box
- Seamless Upholstery
- Two Position Armrests
- 300 LED Light
- 300-Watt Power Supply
- 332 Delivery System with Standard Touchpad
- Doctor Syringe
- Three 4-Hole Vinyl Tubing
- Standard Tray Holder
- Self Contained 2-Liter Water Bottle
- Assistant's Package with HVE, SE and Syringe
- Air and Water QD



### A-dec Headquarters

2601 Crestview Drive  
Newberg, Oregon 97132  
United States  
Tel: 1.800.547.1883 Within USA/CAN  
Tel: 1.503.538.7478 Outside USA/CAN  
Fax: 1.503.538.0276  
a-dec.com/DSO-Solutions  
DSOcustomerservice@a-dec.com


## Choose Your Upgrades



500 LED Light



Cuspidor



Monitor Mount



Lever Foot Control



Deluxe Touchpad



Dr Stool



Cavitron



Camera



Fiber Optic Light Source



A-dec and A-dec logo are registered trademarks in the U.S. Patent and Trademark office. All other products or services mentioned in this document are covered by the trademarks, service marks, or product names as designed by the companies that market those products.

The colors and products shown may vary from the actual materials and are subject to change without notice. Please see your authorized A-dec dealer for accurate color samples and the most current product information.

©2019 A-dec Inc. All rights reserved. 85.0759.00/OP/500/29/Rev C