

A-dec 300 Package

For Dental Service Organizations



 a dec[®]

We've distilled the most essential features into one smart solution. Compact. Sensible. Flexible. A-dec 300 brings you exactly what you need with an added bonus—it's from A-dec, so you know it's built to last.

Features and Benefits

- 60° chair swivel with brake—allows easy rotation and patient positioning.
- Dual-articulating, gliding headrest—supports head and neck, moves with the patient.
- Fully ambidextrous delivery system—change from right-hand to left-hand configuration.
- Hydraulic lift system—A-dec built for robustness and gentle stops and starts during patient positioning.
- Large 2-liter water bottle—requires fewer refills.
- Safety lift arm stop plate—stops immediately on contact for operator safety.
- Contoured floor box—sleek, unobtrusive floor box discreetly houses air and water controls.

Choose Your Upgrades



500 LED Light



Cuspidor



Monitor Mount



Deluxe Touchpad



500 Doctor and Assistant Stools



Ultrasonic Scaler



Clinical products including, fiber optic light source, intraoral camera, motors and handpieces

A-dec 300 DSO Package Includes

- 311 Chair with Contoured Floor Box
- Formed Upholstery
- Two Position Armrests
- 300 LED Light
- 300-Watt Power Supply
- 332 Delivery System with Standard Touchpad
- Doctor autoclavable QD Syringe
- Three 4-Hole Vinyl Tubing
- Standard Tray Holder
- Self Contained 2-Liter Water Bottle
- Assistant's Package with autoclavable HVE, SE and Syringe
- Air and Water QD
- 421 Doctor and 422 Assistant Stools



A-dec Headquarters

2601 Crestview Drive
Newberg, Oregon 97132
United States

Tel: 1.800.547.1883 Within USA/CAN
Tel: 1.503.538.7478 Outside USA/CAN
Fax: 1.503.538.0276
a-dec.com/DSO-Solutions
DSOcustomerservice@a-dec.com

The colors and products shown may vary from the actual materials and are subject to change without notice. Please see your authorized A-dec dealer for accurate color samples and the most current product information.

©2022 A-dec Inc. All rights reserved. 85.0759.00/BB/500/9-22/Rev D