

Performer Package

For Dental Service Organizations



 a dec[®]

Designed to give you flexibility, access and comfort. Performer® answers the need for a lower cost option designed with a rich feature set that supports productivity.

Features and Benefits

- Fully ambidextrous delivery system—change from right-hand to left-hand configuration.
- Contemporary and stylish delivery system.
- User programmable touch pad.
- Hydraulic cantilever system for a smooth quiet operation.
- Robust cast metal construction.
- One-piece floor box.
- Large 2-liter water bottle—requires fewer refills.
- Safety lift arm stop plate—stops immediately on contact for operatory safety.
- Contoured floor box—sleek, unobtrusive floor box discreetly houses air and water controls.
- Standard assistant's instrumentation.

Performer DSO Package Includes

- Chair with One Piece Floor Box
- Seamless Upholstery
- Fixed Armrests
- 300 LED Light
- 300-Watt Power Supply
- Traditional Delivery System with Standard Touchpad
- Doctor Syringe
- Three 4-Hole Vinyl Tubing
- Standard Tray Holder
- Self Contained 2-Liter Water Bottle
- Assistant's Package with HVE, SE and Syringe
- Air and Water QD



A-dec Headquarters

2601 Crestview Drive
Newberg, Oregon 97132
United States

Tel: 1.800.547.1883 Within USA/CAN
Tel: 1.503.538.7478 Outside USA/CAN
Fax: 1.503.538.0276
a-dec.com/DSO
DSOcustomerservice@a-dec.com

Choose Your Upgrades



500 LED Light



Cuspidor



Monitor Mount



Lever Foot Control



Deluxe Touchpad



Dr Stool



Cavitron



Camera



Fiber Optic Light Source



A-dec and A-dec logo are registered trademarks in the U.S. Patent and Trademark office. All other products or services mentioned in this document are covered by the trademarks, service marks, or product names as designed by the companies that market those products.

The colors and products shown may vary from the actual materials and are subject to change without notice. Please see your authorized A-dec dealer for accurate color samples and the most current product information.