



A-dec® 12 o'clock Treatment Packages

Fleet Options*

Tucked away at the head of the chair, the A-dec Inspire 391 Fleet Treatment Console and A-dec 341 Fleet 12 o'clock Duo delivery system create an uncluttered, inviting patient experience. The delivery's dual-user design includes an easy-to-position worksurface and ample housing for up to five integrated clinical devices.



A-dec Inspire 391 Fleet Treatment Console:

- Choice of pass-through or non-pass-through upper storage to align with your treatment room design
- Solid surface countertop in Designer White
- Four-drawer module to maximize in-treatment room organization and storage
- Side-opening access to a single ideally positioned swing-out shelf for the dental assistant
- Quiet soft-close drawers and doors
- Task lighting matches color temperature of the A-dec LED dental light
- USA and Canada UL-approved wiring with hospital-grade GFCI tamper resistant outlets
- Single duplex outlet and USB 2.0
- Glove dispenser
- 8 glass shelves

A-dec 341 Fleet Delivery:

- Standard touchpad supports manual control of chair, light, and cuspidor functions (if equipped), and four programmable chair position presets
- Five adjustable instrument holders
- Three- or four-position handpiece control block with individually adjustable water coolant
- Choose from four easy-to-order configurations
- Control center houses up to five integrated clinical devices
- 2-liter self-contained water system with quick-disconnect (QD) water bottle
- Solids collector, dual high-volume evacuator (HVE), QD syringe, and saliva ejector
- Solid surface 22.5" (571.5 mm) diameter round worksurface in Designer White
- Independent height-adjustment for doctor's instrumentation and worksurface
- L/R delivery positioning supports left- or right-handed operators
- Optional A-dec+ Gateway enables A-dec+ software connectivity

A-dec 341ft Treatment Package Configurations

#1 #2 #3 #4
(77.2270.00) (77.2271.00) (77.2272.00) (77.2273.00)

		#1	#2	#3	#4
Doctor Side*	Radius Cabinet Mount (Inspire 591, 391, Preference Collection 5562,5662,5580)	•	•	•	•
	No Tray Holder Assy	•			
	Standard Tray Holder Assy		•	•	•
	Disc Foot Control, Wet-Dry (only)	•	•	•	•
	Power Supply Connection Kit (341ft to 311/411 chair power supply)	•	•	•	•
	Standard Touch Pad (membrane)	•	•	•	•
	Quad Volt Intraoral Light Source		•	•	•
	Three Position Control Block	•	•		
	Four Position Control Block			•	•
	QD Autoclaveable 3-way Syringe	•	•	•	•
	Pos. 1: No Tubing	•	•	•	•
	Pos. 2: 4-hole Handpiece Tubing	•			
	Pos. 2: 6-pin Handpiece Tubing		•	•	•
	Pos. 3: 4-hole Handpiece Tubing	•	•	•	•
	Pos. 4: 4-hole Handpiece Tubing	•	•	•	•
Integrated 15' USB 2.0 Cable	•	•	•	•	
Assistant Side	Solid Surface 22.5" diameter Round Worksurface: Solid Designer White	•	•	•	•
	Assistant Standard Touchpad				•
	QD Autoclaveable 3-way Syringe	•	•	•	•
	Standard 11mm HVE	•	•	•	•
	Additional 11mm HVE and Holder	•	•	•	•
	Saliva Ejector and Holder	•	•	•	•
	2L Self-Contained Water Bottle	•	•	•	•
**A-dec+ Compatible	•	•	•	•	

*Configured for right-handed operator

**Requires optional A-dec+ Gateway (#43.0531.00)

