



A-dec® 300 Pro 12 o'clock Duo Delivery System

Tucked away at the head of the chair, the A-dec 300 Pro 12 o'clock Duo delivery system creates an uncluttered, inviting patient experience. You'll appreciate the intuitive user interface, easy-to-position worksurface, and ample housing for up to five integrated clinical devices.



341pro shown with CP5i and optional accessories

Key Features:

- Available with CP5 or CP5i interface for doctor; standard touchpad for assistant
- Optional A-dec+ Gateway enables A-dec+ software connectivity
- Five adjustable instrument holders
- Three- or four-position handpiece control block with individually adjustable water coolant
- Control center houses up to five integrated clinical devices
- 2-liter self-contained water system with quick-disconnect (QD) water bottle
- Solids collector, high volume evacuator (HVE), QD syringe, and saliva ejector
- 22.5" (571.5 mm) diameter round work surface (solid surface or quartz)
- L/R delivery positioning supports left- or right-handed operators
- Independent height-adjustment for doctor's instrumentation and worksurface

Control Pad 5 benefits:

- Easy-to-clean glass screen
- Manual chair, light, and cuspidor controls (if equipped)
- Dynamic chair presets
- Four adjustable timers
- Optional scaler integration

Control Pad 5i benefits:

- All benefits of CP5, plus...
- Water coolant vaporizer
- Timed flush
- Air and water coolant on/off/linked
- Optional electric / endodontic integration

MAKE IT A CONNECTED EXPERIENCE!

Add the optional A-dec+ Gateway to access the A-dec+ software platform, a digital foundation for integrating the latest clinical products and future technologies. Explore how A-dec+ can help your practice evolve.

