

A-dec® Choices  
*Which solution is right for you?*

Choices

a<sup>®</sup>dec

# There's an A-dec for everyone.

A-dec offers creative, reliable solutions for practices of every size and specialty. Including yours. No matter which dental chair and delivery system you choose, you'll have peace of mind, knowing that they're designed to work the way you do—and to work without fail. Engineered with legendary quality and reliability, A-dec equipment is tested for 20 years of life, and backed with superior service and support. So you can practice with confidence.



## A-dec 500

*Premier comfort, performance, and connectivity*

Bring intelligent simplicity and infallible peace of mind to your operatory, and superior comfort to your patients.

### A-dec 500 Chair

- Exquisitely smooth and quiet, cradles patients with whole body support
- Impressive range of motion, including our lowest base-down position
- Ultra-thin backrest flexes to create ample legroom under the chair
- Ergonomically designed for optimal access and healthy, natural posture

### A-dec 500 Pro Delivery

- Six adjustable control head positions for customizable placement
- DS7 touchscreen offers a simple, dynamic interface for intuitive, connected control
- Multi-user settings can be individually programmed and personalized
- Enabled with A-dec+ software, to easily add more services and features as technology evolves



## A-dec 400

*Contemporary luxury with proven reliability*

Sleek design, patient comfort, and uninhibited dental team access—all combined with reliable performance.

### A-dec 400 Chair

- Distinctive contemporary style with contoured patient comfort
- Smooth, efficient positioning, plus 30° rotation on each side of center
- Ultra-thin backrest flexes to create ample legroom under the chair
- Ergonomically designed for excellent access and healthy, natural posture

### Delivery Options

- Pairs with A-dec 300 Pro delivery
- Choose from Continental® or Traditional deliveries
- Radius® option rotates easily around the chair to support both left- and right-handed dentistry



## A-dec 300

*A small footprint with full-size performance*

Accessible, flexible, and evolved to bring you rich feature sets and the option of digital connectivity.

### A-dec 300 Chair

- Compact, swiveling design fits a variety of operatory spaces, with minimal maintenance
- Contoured, pressure-mapped surfaces support patients in comfort
- Ultra-thin backrest flexes to create ample legroom under the chair
- Ergonomically designed for excellent access and healthy, natural posture

### A-dec 300 Pro Delivery

- Open platform design is easy to configure, customize, and upgrade in the future
- Pivoting control head reduces hand strain and fatigue
- Intuitive, responsive user interface shows just what you need to know in the moment
- Add the optional A-dec+ Gateway to access the A-dec+ software platform

Chair	A-dec 500	A-dec 400	A-dec 300
Ultra-soft start/stop hydraulic drive	•		
Hybrid drive (hydraulic base/electric tilt)		•	•
Swivel (30° each side of center)	•	•	(US/Canada only)
Virtual pivot for added patient comfort	•	•	•
Thin backrest 1.25" (32mm)	•	•	•
Dual-articulating gliding headrest	Lever release	Locking knob or lever release	Locking knob
Multi-position armrest settings	2	2	2
Chair at lowest point	13.5" (343mm)	13.75" (349mm)	13.75" (349mm)
Chair at highest point	31.5" (800mm)*	31.5" (800mm)	31.5" (800mm)
Integrated 300 watt power supply	•	•	•
Load rating	500 lb. (227 kg)	400 lb. (181 kg)	400 lb. (181 kg)
Choice of sewn or formed upholstery	•	•	•

Delivery System	A-dec 500 Pro	A-dec 300 Pro
Chairs	A-dec 500	A-dec 500, A-dec 400, A-dec 300
Mounting options	Radius and 12 o'clock	Radius, side, and 12 o'clock
Control pad	DS7	CP5, CP5i
A-dec+ connectivity	•	Optional (requires A-dec+ Gateway)
Four-position control block	•	Optional
Three-position control block	—	•
Holder positions	6	5
Clinical device integration	5	2
USB pass-through port	•	(on deliveries with CP5i)
Support center	A-dec 400	A-dec 400, A-dec 300
Self-contained water system with quick-disconnect 2-liter bottle	•	•
Accommodates amalgam separator	•	•
Cuspidor swivel	•	A-dec 400 only
Programmable timing for cup and bowl	•	•
Dental light	A-dec 500 LED	A-dec 500 LED, A-dec 300 LED

\*Using maximum upper limit. The factory set upper limit is 29" (737mm).



## More than a warranty. It's a commitment.

At A-dec, our focus is on making dentistry better—which includes superior customer service along with superior dental products. Your A-dec equipment is backed by a five-year warranty, and a lifetime of customer support. We'll be here today, tomorrow, and well into the future. That's our commitment to you.



### **A-dec Headquarters**

2601 Crestview Drive  
Newberg, Oregon 97132 USA  
**P:** 1.800.547.1883 within USA/CAN  
**P:** +1.503.538.7478 outside USA/CAN  
**F:** +1.503.538.0276  
**a-dec.com**

### **Other Locations:**

A-dec Nashville Experience Center  
A-dec Australia  
A-dec China  
A-dec United Kingdom

Certain touchpad symbols are proprietary to A-dec Inc. Any use of these symbols, in whole or in part, without the express written consent of A-dec Inc. is strictly prohibited.

The colors and products shown may vary from the actual materials and are subject to change without notice. Please see your authorized A-dec dealer for accurate color samples and the most current product information.

©2023 A-dec Inc. All rights reserved.  
85.0394.00/BB/1.5M/8-23/RevM