

Maximize your vacuum performance with A-dec®.

A-dec Evacuation System Cleaner

A vacuum pump is at the core of a practice. Over time, organic and inorganic materials may build up in the vacuum lines and clog traps, impacting the performance of your vacuum system.

Keep your pump running at peak performance with A-dec Evacuation System Cleaner to avoid costly repairs and downtime. A-dec's unique, multi-enzyme solution is specifically formulated to dissolve and remove organic and inorganic debris that may build up over time.

Applying the solution is simple. Mix one ounce (30 mL) with one liter of water and dispense the solution using the A-dec ICV® vacuum line cleaning system or an atomizer.

WHY A-DEC?

- Multi-enzyme solution
- Non-foaming, pH neutral
- Compatible with amalgam separators
- Removes organic and inorganic debris
- Fresh scent eliminates odors
- Compatible with A-dec vacuum pump and ICV vacuum line cleaning system



Product Description	Contents	A-dec Part Number
33 fl oz (976 mL) Liquid Concentrate	12 Bottles	91.0060.00
33 fl oz (976 mL) Liquid Concentrate	Single Bottle	91.0060.01
128 fl oz (3.78 L) Liquid Concentrate	4 Bottles with Pumps	91.0061.00
128 fl oz (3.78 L) Liquid Concentrate	Single Bottle with Pump	91.0061.01

To order A-dec Evacuation System Cleaner, contact your local dealer or request a quote at a-dec.com/linecleaner