



Increase aerosol capture without replacing your equipment.



A-dec® Large Bore HVE

The better you control aerosols at the treatment site, the lower the risk of viral exposure. The A-dec® Large Bore HVE removes 22-45% more air from the oral cavity than our standard HVE. It easily integrates into your existing A-dec equipment with no additional tubing, connections, or operatory space needed. It's one small change that can make a big difference for your peace of mind.



- Removes 22-45% more air volume than our standard HVE
- Adaptable to corresponding ICV system
- Integrates into current A-dec equipment
- Functional for single-operator use
- Fast and easy to install
- Tips are fully autoclavable

Product Description	Contents	Equipment Compatibility	A-dec Part Number
Large Bore HVE Kit	15mm HVE with 63" tubing & holder	352, 353, 541 & 545	11.1370.00
Large Bore HVE Kit	15mm HVE with 63" tubing, holder & 20 tips	352, 353, 541 & 545	11.1371.00
Large Bore HVE Kit	15mm HVE with 52" tubing & holder	351 & 551	11.1372.00
Large Bore HVE Kit	15mm HVE with 52" tubing, holder & 20 tips	351 & 551	11.1373.00
Dürr Universal Cannula III	20 tips	N/A	11.1374.00

To learn more about the A-dec Large Bore HVE, visit a-dec.com/360 or contact your local authorized dealer.