

A-dec PACI

Portable Dental System



 a dec[®]

Durable. Simple. Versatile.

The PAC I portable dental unit is field-proven to go wherever you need to be.

- Manual control for two handpieces
- Wet/dry foot control
- Non-retracting handpiece water coolant
- Water and vacuum systems that quickly clip on and off
- Self-contained water system
- Storage within the control head for handpiece tubing, syringe, clean water system and accessory air supply line
- Complete repair kit for quick servicing anywhere
- Unit top and housing surfaces that are smooth for easy clean-up
- Side-mounted controls that are out of the aerosol zone
- Smooth, cleanable handpiece tubing
- Autoclavable syringe with smooth tubing and a streamlined design for fast clean-up
- Instrument holders that are easy to clean and are designed for barrier protection
- Handpiece Flush System to flush handpieces quickly and simply



“With the Pac I, we get a dependable workhorse. With A-dec, we have a partner that supports us with whatever we need.”

Bruce Bergstrom
America's Dentists Care Foundation Missions of Mercy

Field Unit

The PAC I Field Unit control housing mounts onto an adjustable-height tripod. The Field Unit is specifically designed for access to remote locations and includes its own carrying case.* The unit, optional case, and standard accessories weigh a combined 35 pounds—a light load for easy transport. And if you need to make an adjustment or repair out in the field, a special tool kit is included.

A-dec's Steril-Vac assistant's instrumentation's one-quart capacity waste container is standard with the Field Unit.

Institutional Unit

The Institutional Unit, with its adjustable-height mobile U-frame, is intended for hospitals and other medical care facilities.*

The Institutional Unit includes all the features of the PAC I control and—like the Field Unit—comes standard with the Steril-Vac assistant's instrumentation.

Self-Contained Unit

The PAC I Self-Contained Unit includes its own air supply system. Integrated with the fixed-height mobile frame is a 1/2-hp, 120 or 240 VAC air compressor with a four-liter reservoir. A quiet cooling fan keeps the compressor running efficiently.

The specially designed air saliva ejector is compact, requires less air, and is easily removed for cleaning.

[*] Requires a compressed air source.



The PAC I is the most portable operator system in the field, and the handpiece delivery control housing is available in three versions: Field, Institutional, and Self-Contained.

A-dec Porta-Chair 3460

Designed for high-performance portability. Produced with quality alloy construction and durable upholstery, this portable chair is ultra-light, folds quickly, and—if needed—fits inside an optional carrying case for travel and storage.

- Height adjustment crank (mounts for left or right)
- Adjustable chair heights: seat level 200 mm-710 mm (8" to 28") headrest level 200 mm-1420 mm (8" to 56")
- Thin back (tilts from 75° upright to full horizontal)
- Arm slings
- Adjustable U-shape headrest
- Polyurethane foam and heavy weight upholstery
- Quick disconnect toeboard for compact packing
- Optional zipper carrying bag

To learn more about A-dec portable dental solutions, contact your authorized A-dec dealer.



A-dec Headquarters

2601 Crestview Drive
Newberg, Oregon 97132 USA
Tel: 1.800.547.1883 within USA/CAN
Tel: 1.503.538.7478 outside USA/CAN
Fax: 1.503.538.0276
a-dec.com

A-dec Australia

Tel: 1.800.225.010 within AUS
Tel: +61.(0)2.8332.4000 outside AUS

A-dec China

Tel: +1.503.538.7478

A-dec United Kingdom

Tel: 0800.233285 within UK
Tel: +44.(0).24.7635.0901 outside UK

A-dec, the A-dec logo, A-dec 500, A-dec 300, Cascade, Cascade Master Series, Century Plus, Continental, Decade, ICX, ICV, Performer, Preference, Preference Collection, Preference ICC, and Radius are trademarks of A-dec Inc. and are registered in the United States and other countries. A-dec 400, A-dec 200, Preference Slimline, and reliablecreativesolutions are also trademarks of A-dec Inc. None of the trademarks or trade names in this document may be reproduced, copied, or manipulated in any manner without the express, written approval of the trademark owner.

The colors and products shown may vary from the actual materials and are subject to change without notice. Please see your authorized A-dec dealer for accurate color samples and the most current product information.

©2016 A-dec Inc. All rights reserved.
85.6008.00/OP/1.5M/5-16/Rev C