

A-dec®

Specialty Chair

Get all the long-lasting performance and reliability of A-dec, in a chair made especially for your practice. Together with the A-dec Specialty Workstation, this fixed-base chair ensures a smooth, efficient treatment flow—and a comfortable experience for both you and your patients. Choose from two distinctive upholstery styles, and customize the color to perfectly compliment your practice design.

317 Specialty Chair

Fixed-base chair for orthodontic, pediatric, and specialty practices

Standard Features:

- Compact, space-saving design
- Dual-articulating gliding headrest for optimal oral cavity positioning
- Virtual pivot moves with patient, eliminating need to readjust chair after reclining or returning
- Pressure mapped, contoured surfaces offer lasting comfort and support
- Armless design for easy patient entry and exit
- Thin, flexible backrest lets you roll in close to the patient
- Chair recline switch with automatic back up
- Fixable height setting from 19.5" to 16.5" (495mm to 419mm)
- Supports up to 350 lbs (159kg)*
- Backed with a 5-year A-dec warranty

Options:

- Upholstery styles - sewn or formed (shown)
- Mounting for A-dec LED lights (371L or 571L)

To learn more about the A-dec Specialty Chair, contact your A-dec partner or local authorized dealer.

a-dec.com/partner

*Load rating based on UL requirement of 4X (1400 lb/635 kg) overload static

©2024 A-dec Inc. All rights reserved. 85.1069.00/E-only/8-24/Rev A

