

## Specialty Workstation

Ideal for orthodontics and pediatrics, this versatile design offers easy access and efficient organization for technology and materials, while the articulating arm keeps handpieces comfortably within reach. Customize the surfaces to suit your style—and enjoy the high quality and reliable performance of an A-dec® solution made for your practice.

### A-dec Inspire® 397 Cabinet

#### Standard Features:

- Stationary base
- 28" model measurements:  
28"W x 32 5/8"D x 32"H (right hand only)
- 24" model measurements:  
24"W x 32 5/8"D x 32"H (right hand only)
- 1 cm solid surface countertop
- Five full-extension 5" drawers
- Quiet soft-close drawers
- Computer storage with fan
- Subbase with easy access to utilities
- 341sp delivery system accommodation
- USA and Canada UL-approved wiring with hospital grade tamper-resistant outlets
- Duplex convenience outlet
- Quad outlet in subbase

#### Optional Upgrades:

- 2 cm quartz countertop
- A-dec LED dental light and monitor mount accommodation
- ICV accommodation
- Glove dispensers (1 or 2)
- Waste drop
- Wipes holder

### A-dec 341sp Delivery System

#### Standard Features:

- Standard touchpad
- 5 position holder  
Position 1: 4-hole  
Position 2: 4-hole  
Position 3: Syringe  
Position 4: SE  
Position 5: HVE
- Wet/dry foot control
- Air and water QDs
- Chair power supply connector kit

### Monitor & Light Options

#### Light:

- A-dec 378L Specialty Dental Light\*
- A-dec 578L Specialty Dental Light\*

#### Monitor Mount:

- 381 monitor mount
- 482 monitor mount
- Optional monitor mount handle

*\*includes power supply and lower post*

Workstation shown with optional upgrades.

©2023 A-dec Inc. All rights reserved.

85.1056.00/E-only/1-24/Rev B3

To learn more about the Specialty Workstation,  
contact your local authorized dealer.