

# Air Motor Sterilization? Now it's a snap.



## Quick Disconnect Air Motors from A-dec®

Dentistry has more safety guidelines and requirements than ever before. Here's a smart, simple solution for best-practice sterilization. **W&H® Quick Disconnect Air Motors** detach from handpiece couplers with minimal effort, to allow thorough cleaning and sterilizing between patients. Fast. Easy. Efficient. For compliance without compromise.

- Quick disconnect motor/coupler design eliminates cumbersome threading
- Sleek inline design is easier to use than large, bulky motors
- Slim and lightweight to reduce hand and wrist fatigue
- High torque ensures faster, more efficient cutting
- Smooth, quiet operation puts patients at ease
- Compatible with popular Midwest® and StarDental® brands\*
- Backed with 2-year A-dec warranty
- Requires RQ-04 coupler for non-lighted handpiece



Product Name	Handpiece Compatibility*	Max Speed (rpm)	Sterilizable	A-dec Part Number
AM-20RQ	A-dec®   W&H®	20,000	Yes	0.30400000
AM-20MWRQ	Midwest®	20,000	Yes	0.30402000
AM-5MWRQ	Midwest®	5,000	Yes	0.30404000
AM-20SDRQ	StarDental®	20,000	Yes	0.30401000
AM-5SDRQ	StarDental®	5,000	Yes	0.30403000

To learn more about Quick Disconnect Air Motors, contact your local authorized dealer.



Not available in all countries. Check for local availability.

\* Midwest® and StarDental® are trademarks of Midwest Dental and DentalEZ respectively. Midwest Dental and DentalEZ are not affiliated with W&H or A-dec.  
©2020 A-dec, Inc. All rights reserved. 85.0889.00/6-20/Rev A

